

Vivant

Range Modifications

Release date: August 1st 2001

WHEN QUALITY MATTERS

The Vivant product has undergone some design changes in order to enhance build quality. Customer valued features & aesthetics have not been changed.

Body

Vivant body construction has changed – becoming a one piece folded steel design.

Fully integral emergency gear has moved from the centre compartment, with painted metal cover, to a back access removable gear tray mounted behind the profiler.

The body is still compatible with the following ceilings:

- 15 & 25mm Exposed Tee.
- 15mm Tegular & Decorative Tee.
- Spring Tee.
- Cut aperture.




Profilers

The luminance profiler system, projecting light onto the ceiling, has been redesigned – improving the ceiling illumination.

- The opal plastic extruded upper diffuser, fixed behind/above the profilers, has been replaced by folded steel, painted to match the body.

Vivant is still available with Lift, Matrix & Blend profilers.



Louvre

Slight modifications have been made to the louvre.

- The louvre is now offered with 13 wedge blades for Retail and 15 wedge blades for Cat2 options.
- The louvre centre wedge now covers the control gear - eliminating the need for a separate cover.
- Aluminium louvre material has been upgraded from 610G3 to MIRO6.
- Louvre now comes fitted with dog-tag safety cords.

